RGP CONTACT LENS INSTRUCTIONS FOR USE

THE END USER OF THIS DEVICE MUST READ THESE INSTRUCTIONS FOR USE

CSD.81 (ISSUE 12 - FEBRUARY 2021)
(ORIGINAL ISSUE - JULY 1996)

EMERGO
EUROPE

RGP PRODUCTS COVERED BY THIS INSTRUCTIONS FOR USE

4F

APEX

AQUALINE MF

CUSTOM

INTEGROLENS

METRO2

NISSEL EIGHTEEN-20

NISSEL SIXTEEN-50

NISSEL KII RIGID

FOR MORE INFORMATION ABOUT:

**OUR RGP
PRODUCTS, SCAN**

**OUR INSTRUCTIONS
FOR USE, SCAN**

**OUR CERTIFICATION,
SCAN**

YOUR RGP LABEL

YOUR LENS PARAMETERS

YOUR AUTOMATED SPECIALISED LOT NUMBERS

R

L

LOT

LOT

ORDER

RE-SANITISE AFTER 28 DAYS.

MD

! 

CANTOR & NISSEL LTD
(MARKET PLACE,
BRACKLEY,
NORTHAMPTONSHIRE,
NN13 7NN, UK)

SINGLE
PATIENT
USE ONLY

CE

1639

EC REP

EMERGO EUROPE
(PRINSESSEGRACHT 20,
2514 AP THE HAGUE,
THE NETHERLANDS)

YOUR SANITISED DATE AND ORDER NUMBER

THE MANUFACTURER

CAUTION, CONSULT ACCOMPANYING DOCUMENTATION

MEDICAL DEVICE

EC REP [OUR EC REPRESENTATIVE, EMERGO EUROPE]

RGP CONTACT LENS - INSTRUCTIONS FOR USE

Lenses should not be worn unless prescribed by a registered eye care professional. The end user of the device **MUST** read these instructions **BEFORE** use – important information is inside.

Description

All RGP contact lenses manufactured by Cantor & Nissel Ltd are covered under this instructions for use.

Intended Use

Correction of vision for conditions such as Myopia, Hyperopia, Presbyopia, Astigmatism & Keratoconus.

Indications

Spherical lens are indicated for myopia, hyperopia. Aspherical bifocal lens is indicated for myopia, hyperopia, and presbyopia. Toric lenses are indicated for myopia, hyperopia, astigmatism. Certain lenses are indicated for use with Keratoconus.

Instructions For Use

These lenses are daily wear rigid contact lenses. The lenses wear time will be suggested by your eye care professional. Your eye care professional will recommend lens wearing, cleaning regime, and replacement regime. Do not vary these. If in doubt seek the advice of your eye care professional.

Instructions on Handling Lenses & Contact Lens Solutions

1. Keep your finger nails short and smooth to avoid scratching your eyes or your lenses.
2. Always wash your hands with warm water and soap and dry thoroughly before handling your lenses. Lenses may become contaminated if you do not wash your hands properly.
3. Always handle your lenses somewhere that is well-lit, clean and sterile.
4. Your contact lenses may become damaged if mishandled.
5. The risk of infection is higher without proper care for your contact lenses and cleaning solutions

Checking the lenses and packaging

1. Check which lenses are for the right eye and which lenses are for the left eye. (if applicable)
2. Check if there is anything wrong with the lenses, plastic packaging, or tamper evident seal. If there is anything wrong, for example, if there is any damage or leakage, do not open or use and contact the point of purchase.

Removing the lenses from the plastic packaging

1. Hold the package and use one hand to pull upwards against the tamper evident seal of the packaging.
2. Remove the packaging from the plastic tamper evident sealed bag and carefully pop open the lens packaging.

Checking the lenses

1. Hold the edges of the lens from the side and gently move it up and down to remove the liquid.
2. Put the convex surface on the index finger of your dominant hand and check the lens (Fig. 1)
3. Check if there is anything wrong with the lens, for example, if the lens is deformed, discoloured, contains foreign matter, is damaged, is scratched or is dirty. Do not wear the lens if you notice any abnormalities.

Inserting Lenses

1. Put the lens onto your finger, and with your middle finger of the same hand, pull down your lower eyelid. With your index finger of your other hand, pull up your top eyelid. If you put your finger on your eyelash line, you can open your eye wide.
2. Look in the mirror so that your cornea is located centrally. Slowly bring the lens close to your eye and put it gently in place whilst still looking in the mirror. If you glance upwards when doing this, the lens will not be properly centered on your cornea.
3. When you have placed the lens properly on the cornea, slowly remove your index finger and then release your eyelids, starting with the bottom lid first and then the top lid. If you remove your fingers too quickly or blink too quickly, the lens may change position or fall out.
4. Check in the mirror that the lens is on top of your cornea and covering the opposite eye with your hand, check your sight in both eyes. Follow the same procedure to fit the lens in your other eye.

Figure 1

Figure 2

Figure 3

Removing the lenses - Figure 2 & 3

1. Check the position of the lens in the mirror.
2. Looking straight in the mirror, tilt your head slightly down whilst still looking in the mirror.
3. Using both index fingers place them on both upper and lower lids.
4. In this position, use your fingers and lids to very carefully prize the lens from the surface of your cornea. You can also open your eye wide and use one index finger on the corner of your eye, pull the skin to the side and blink. This should pop the lens from the surface of the cornea utilizing your lids.

Contact Lens & Contact Lens Case Care

1. Since proteins and lipids in your tears will become attached to the surface of your contact lenses as you wear them, it is important to follow the correct lens care regime for washing, rinsing, disinfecting and storing your lenses.
2. Failure to follow the lens care regime can cause lenses to become dirty and may result in a serious eye disorder or eye infection.
3. Use lens care products suitable for RGP contact lenses. Only use the recommended chemical cleaning solution. RGP lenses should not be autoclaved.
4. Use lens care products according to their accompanying booklet or instructions
5. The extent to which lenses become dirty varies from person to person. If your lenses tend to become dirty, try rubbing your lenses with a cleaning solution or using a protein removal liquid. Dispose of the lenses if you are unable to clean them sufficiently.
6. Replace your lens case roughly every one or two months. Prolonged use will result in a proliferation of bacteria due to the build-up of dirt, rendering disinfection inadequate and potentially leading to eye infection.
7. If your contact lenses have not been worn for a length of time you must re-sanitize your contact lenses for the full contact time of the solution used and use a new contact lens storage case.

Contact Lens Storage and Transportation

1. Your contact lenses will remain sanitized (see sanitized date on label) until you break the tamper evident seal. If the lenses are not worn directly after doing this, your lenses must be stored straight away in an appropriate contact lens case with contact lens solution.

2. Contact lenses can be transported in their original plastic packaging. Contact lenses which have been removed from their packaging should be transported in an adequate contact lens case to ensure that they stay sanitized and to prevent the ingress of dirt and bacteria.

Always use fresh, unexpired lens care solutions that are recommended by your eyecare professional.

Contradictions

- Medical Conditions: acute and sub-acute inflammation of the anterior eye, eye infection, uveitis, reduced corneal sensitivity, dry eye or lacrimal apparatus disease to an extent that poses a problem for wearing lenses, eyelid disorders, or allergy to an extent that affects the wearing of lenses, and other diseases which would complicate the use of contact lenses.
- Environmental Conditions: anyone living or working in dry environments. Anyone living or working in an environment where dust or chemicals are likely to enter the eye.
- Personal Qualities: anyone who cannot follow the recommendations of their eye care professional, anyone who cannot use contact lenses properly, anyone who fails to attend regular check-ups, anyone who cannot follow the hygiene routine required to wear contact lenses.
- Examples of uncommon eye disorders that can be caused by RGP contact lenses: corneal oedema, corneal ulcer, keratitis (including infectious keratitis & acanthamoeba keratitis), conjunctivitis, corneal abrasion and corneal neovascularization

Lenses may be damaged or become damaged. If you have inadvertently worn a lens which is damaged or otherwise faulty, or if a lens become damaged while you are wearing it, remove the lens immediately and consult your eye care professional.

Precautions

Be sure to read this booklet carefully before using your contact lenses, and if you are in any doubt about the terms or information contained herein, please consult your eye care professional for clarification before using the medical device. Please retain this document for future reference.

- Do not wear contact lenses that are damaged.
- Do not wear contact lenses for longer than the recommended wearing period.
- Do use contact lenses properly in accordance with recommended methods for handling contact lenses.
- Always attend regular contact lens check-ups.
- If you experience any problems at all, consult your eye care professional immediately.
- Always insert lenses promptly after opening, as lenses can become easily contaminated. If you do not use them straight away, store them according to the instructions of the lens care products used.
- Avoid rubbing your eyes while wearing lenses, as this may result in scratches to the surface of your cornea or lens displacement.
- Do not mix up right and left lenses, as doing so may result in impairment of vision and discomfort.
- Make sure you have appropriate glasses to use in case you are unable to wear your contact lenses due to a condition of your eyes or for other reasons, or after you have worn your lenses longer than the recommended wearing period stated by the eye care professional. Continuing to wear your contact lenses regardless may cause eye disorders or other problems.
- Follow the instructions of your eye care professional if using eye drops while wearing contact lenses.
- Do not allow cosmetics, pharmaceuticals or other substances to come into contact with your lenses. If cosmetics (eyeliner, cleansing products, etc.) contaminate your lenses this could result in discoloration, irritation and/or deterioration of the quality of the contact lens.
- When your eyes and contact lenses are exposed to strong wind (such as riding a bicycle, skiing, or direct wind etc.) always use some kind of protection against the wind such as sunglasses or goggles. Failure to do so may result in uncomfortably dry lenses and/or potential displacement of lens.

- If solvents such as household detergents or cleaning products get in your eyes, remove and dispose of your lenses immediately, and consult your eye care professional.
- Never wear someone else's contact lenses or give your lenses to someone else to wear, as this could result in possible infection and serious eye disorders.
- Users who find it difficult to insert lenses, such as the elderly, young children or severely vision impaired should first consult an eye care professional for guidance.
- Some pregnant women may find that lenses (which were once comfortable) may have become uncomfortable and experience a temporary deterioration of vision. Guidance on further lens wear should be consulted with an eye care professional.
- People who suffer from allergies may be more prone to diseases of the eye and skin than other wearers of contact lenses. Guidance should be received from an eye care professional before wearing of contact lenses.
- Use of contact lenses may be affected by illness, pregnancy, childbirth, use of medications, physical conditions, eye conditions and living environment. The conditions in which you wear contact lenses may also change after you have started wearing contact lenses.

If you experience any pain or continued discomfort, discontinue wear and consult an eyecare professional immediately.

Warnings

- To wear your contact lenses safely, be sure to read this booklet carefully and to follow the instructions of your eye care professional and the proper methods for handling contact lenses.
- Wearing contact lenses and the direct exposure of contact lenses to non-sterile water sources (e.g. tap water, bathing, showering, swimming, participating in waters sports etc.) severely increases the risk of microbial infection and permanent eye damage.
- The use of non-sterile water (tap water, etc) in the handling of contact lenses and contact lens cases severely increases the risk of microbial infection and permanent eye damage.
- Do not change contact lens care systems, lens types or parameters without first consulting an eyecare professional. Improper handling and use of contact lens solutions can result in serious infection. Follow all instructions given by the manufacturer of the contact lens solution, making sure to disinfect lenses for the full contact time stated by the solution used.

Potential Adverse Effects

Until you are used to your lenses, you may experience symptoms such as the following (below). Symptoms vary from person to person but will gradually diminish as you become used to your lenses. However, if you experience acute symptoms or if symptoms persist, stop wearing your lenses and consult your eye care professional.

- Slight discomfort
- Unstable vision
- Watery eyes
- Itchy eyes

Disposal of Contact Lenses

Contact lenses can be disposed of in standard waste. Contact lenses should NOT be flushed down the toilet, or disposed of down the sink. The preferred method of disposing of contact lenses is through a dedicated contact lens recycling program.

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